


# DigiRad

digital radiography

## Capture & View EQUINE


Portable Optimized Imaging Software & Workstation



Software meets Normative and veterinary codes to achieve high **Performance**




Commitment to supporting equine veterinarian with **Innovative tools** for diagnostic and high-end imaging



**Configurable** and orderly pre-purchase with automatism in the preparation process and saved with rectifications



Acquisition from FPDs applying the last generation of algorithms with the best initial image in **< 5 secs**




Best image quality with independence of x-ray source used, to achieve high **Productivity** for any equine workflow



Just a single screen for acquisition process including FPD, X-ray and APR control tools

All images of the study on screen including original and processed

Fully internationally approved APR including zone, subzone and view



Worklist

Complete management with emergency

Capture Module

Main Equine FPD's integrated

Image Tools

Window / Level, Pan, Zoom, Flip, Rotate, Invert, Aut LUT, annotations, angle, measurements, online adjustable magnify glass

Database

SQL equine DB

Export Patients

Export of selected images of a patient to the required file format/s

Images output

DICOM 3, TIFF, JPEG (w/DICOM Tags) BMP, Mail, PACS

Pre-Purchase

Automatic unvalidated workflow on laptop. Configurable & unlimited Reviewable in process time. Autosave




Image Processing

Intelligent Image Processing for equine images. Image processing algorithms configurable by user and applicable to specific zone  
Predefined image processing modes with different intensities

Image Area

25 x 30 cm

AED

Auto trigger

Pixel size

127µm

Spatial resolution

4.2 lp/mm

Scintillator

CsI: TI

Number of bits

16 bit; 65536 Grays

LAN

Gigabit Ethernet

WiFi

802.11n WiFi

Batteries

Rechargeable Lithium Ion 4.5 h

Weight

1.68 kg incl. bat

Warranty

2 years



1012WCC

Tablet and Laptop

Intel i5, 3.2 GHz

MS-Windows 10 Prof. 64 bit

500 GB - 1 TB

4 - 8 GB DDR4

Multitouch LED+ HD 1920x1080

Gigabit Ethernet

Dual AC 802.11n

3 years

Computer

Operating system

Storage

RAM


Screen

LAN

WiFi

Warranty

## Accesorios



**DigiRad**  
digital radiography